

ANDOCOR



ARTERIAL CANNULAE



Reinforced, Curved Tip With Flange

	without vent plug	with vent plug	size french	length cm	qty / master carton
without connector	A2016	A2016V	20	30	20
	B2016	B2016V	20	21	20
	A2216	A2216V	22	30	20
	B2216	B2216V	22	21	20
	A2416	A2416V	24	30	20
	B2416	B2416V	24	21	20
with connector	A20161	A20161V	20	30	20
	B20161	B20161V	20	21	20
	A22161	A22161V	22	30	20
	B22161	B22161V	22	21	20
	A24161	A24161V	24	30	20
	B24161	B24161V	24	21	20
with luer connector	A20162	A20162V	20	30	20
	B20162	B20162V	20	21	20
	A22162	A22162V	22	30	20
	B22162	B22162V	22	21	20
	A24162	A24162V	24	30	20
	B24162	B24162V	24	21	20

ARTERIAL CANNULAE



Reinforced, Straight Tip with Ring

	without vent plug	with vent plug	size french	length cm	qty / master carton
without connector	A2013	A2013V	20	30	20
	B2013	B2013V	20	21	20
	A2213	A2213V	22	30	20
	B2213	B2213V	22	21	20
	A2413	A2413V	24	30	20
	B2413	B2413V	24	21	20
with connector	A20131	A20131V	20	30	20
	B20131	B20131V	20	21	20
	A22131	A22131V	22	30	20
	B22131	B22131V	22	21	20
	A24131	A24131V	24	30	20
	B24131	B24131V	24	21	20
with luer connector	A20132	A20132V	20	30	20
	B20132	B20132V	20	21	20
	A22132	A22132V	22	30	20
	B22132	B22132V	22	21	20
	A24132	A24132V	24	30	20
	B24132	B24132V	24	21	20

4
version 03 - dec 09

ARTERIAL CANNULAE



Non- Reinforced, Curved Tip with Flange

	without vent plug	with vent plug	size french	length cm	qty / master carton
without connector	A2006		20	22	20
	A2206		22	22	20
	A2406		24	22	20
with connector	A20061	A20061V	20	22	20
	A22061	A22061V	22	22	20
	A24061	A24061V	24	22	20
with luer connector	A20062	A20062V	20	22	20
	A22062	A22062V	22	22	20
	A24062	A24062V	24	22	20

5
version 03 - dec 09

ARTERIAL CANNULAE



Non- Reinforced, Straight Tip with Ring

	without vent plug	with vent plug	size french	length cm	qty / master carton
without connector	A2003		20	22	20
	A2203		22	22	20
	A2403		24	22	20
with connector	A20031	A20031V	20	22	20
	A22031	A22031V	22	22	20
	A24031	A24031V	24	22	20
with luer connector	A20032	A20032V	20	22	20
	A22032	A22032V	22	22	20
	A24032	A24032V	24	22	20

AORTIC CATHETER



Non- Reinforced, Straight Tip with Ring

	without vent plug	with vent plug	size french	length cm	qty / master carton
without connector	AC2013		20	21	20
	AC2213		22	21	20
	AC2413		24	21	20
	AC2613		26	21	20
	AC3013		30	21	20
with connector	AC201341	AC201341V	20	21	20
	AC221341	AC221341V	22	21	20
	AC241341	AC241341V	24	21	20
	AC261341	AC261341V	26	21	20
	AC301341	AC301341V	30	21	20
with luer connector	AC201342	AC201342V	20	21	20
	AC221342	AC221342V	22	21	20
	AC241342	AC241342V	24	21	20
	AC261342	AC261342V	26	21	20
	AC301342	AC301342V	30	21	20

VENOUS CANNULAE



Straight, Wire Reinforced

rigid tip	size french	length cm	qty / master carton
V2018	20	39	20
V2618	26	39	20
V3218	32	39	20
V3418	34	39	20
V3618	36	39	20
V4018	40	39	20

Right Angled, Wire Reinforced

V2019	20	39	20
V2619	26	39	20
V3219	32	39	20
V3419	34	39	20
V3619	36	39	20
V4019	40	39	20

VENOUS CANNULAE



Non-Reinforced

straight tip	connection	size french	length cm	qty / master carton
V3208	3/8"	32	50	20

TWO STAGE VENOUS CANNULAE



With Rigid Tip
Proximal And Distal Wire Reinforced

	rigid tip	flat body	size french	length cm	qty / master carton
without connector	V321401	V321401FB	32/40	40	20
	V341401	V341401FB	34/40	40	20
	V361511	V361511FB	36/51	40	20
with connector	V3214011	V3214011FB	32/40	40	20
	V3414011	V3414011FB	34/40	40	20
	V3615111	V3615111FB	36/51	40	20
with luer connector	V3214012	V3214012FB	32/40	40	20
	V3414012	V3414012FB	34/40	40	20
	V3615112	V3615112FB	36/51	40	20

TWO STAGE VENOUS CANNULAE



With Rigid Tip
Proximal Reinforced

	rigid tip	flat body	size french	length cm	qty / master carton
without connector	V360511	V360511FB	36/51	40	20
with connector	V3605111	V3605111FB	36/51	40	20

Distal Reinforced

without connector	V361510		36/51	40	20
with connector	V3615101		36/51	40	20

Non Reinforced

without connector	V360510		36/51	40	20
with connector	V3605101		36/51	40	20
with luer connector	V3605102		36/51	40	20

FLEX-LINE VENOUS CANNULAE



Reinforced with Lighthouse Tip

straight tip	right angled tip	size french	length cm	qty / master carton
01V321L8	01V321L9	32	39	20
01V361L8	01V361L9	36	39	20

FLEX-LINE TWO STAGE VENOUS CANNULAE



Proximal and Distal Reinforced with Lighthouse Tip

	standard body	flat body	size french	length cm	qty / master carton
without connector	01V32L40	01V32L40FB	32/40Fr	40	20
	01V34L46	01V34L46FB	34/46Fr	40	20
	01V36L46	01V36L46FB	36/46Fr	40	20
	01V36L51	01V36L51FB	36/51Fr	40	20
with connector	01V32L401	01V32L401FB	32/40Fr	40	20
	01V34L461	01V34L461FB	34/46Fr	40	20
	01V36L461	01V36L461FB	36/46Fr	40	20
	01V36L511	01V36L511FB	36/51Fr	40	20
with luer lock	01V32L402	01V32L402FB	32/40Fr	40	20
	01V34L462	01V34L462FB	34/46Fr	40	20
	01V36L462	01V36L462FB	36/46Fr	40	20
	01V36L512	01V36L512FB	36/51Fr	40	20

FLEX-LINE TWO STAGE VENOUS CANNULAE



Proximal and Distal Reinforced with Bullet Tip

	standard body	flat body	size french	length cm	qty / master carton
without connector	01V34140	01V34140FB	34/40Fr	40	20
	01V34146	01V34146FB	34/46Fr	40	20
	01V36146	01V36146FB	36/46Fr	40	20
	01V36151	01V36151FB	36/51Fr	40	20
with connector	01V341401	01V341401FB	34/40Fr	40	20
	01V341461	01V341461FB	34/46Fr	40	20
	01V361461	01V361461FB	36/46Fr	40	20
	01V361511	01V361511FB	36/51Fr	40	20
with luer lock	01V341402	01V341402FB	34/40Fr	40	20
	01V341462	01V341462FB	34/46Fr	40	20
	01V361462	01V361462FB	36/46Fr	40	20
	01V361512	01V361512FB	36/51Fr	40	20

CARDIOPLEGIA CANNULAE AND CATHETERS



Aortic Root Cannulae, Non Vented

catalog #	size needle	size french	length cm	qty / master carton
ARC07	18G	7	15	20
ARC04	16G	8	15	20
ARC01	14G	10	15	20

Aortic Root Cannulae, Vented

ARC08	18G	7	15	20
ARC05	16G	8	15	20
ARC02	14G	10	15	20

CARDIOPLEGIA CANNULAE AND CATHETERS



Retrograde Cardioplegia

catalog #	size needle	qty / master carton
RCC01	14G	20



Cardioplegia Needle

CN16	16G	20
------	-----	----

CARDIOPLEGIA ADAPTERS



Cardioplegia Adapters

description	catalog #	qty / master carton
Coronary perfusion Y adapter with male luer locks, female luer lock infusion connector, color-coded clamps	CSY03	20
Multiple perfusion set for simultaneous perfusion of aortic root and three vein grafts with male luer locks, female luer lock infusion connector, color-coded clamps	CSY04 photo not available	20
Multiple perfusion set for simultaneous perfusion of aortic root and three vein grafts with male luer locks and attached vessel cannula, female luer lock infusion connector, color-coded clamps	CSY05 photo not available	20
Multiple perfusion set for simultaneous perfusion of aortic root with vent and three vein grafts with male luer locks, female luer lock infusion connector, color-coded clamps	CSY07 photo not available	20
Aortic perfusion set (with vent), female luer lock, male luer lock infusion connector, color-coded clamps	CSY08	20


CARDIOPLEGIA ADAPTERS



CSY09



CSY10



CSY11



CSY12

Cardioplegia Adapters

description	catalog #	qty / master carton
Cardioplegia set for dual solution administration with female luer locks, male luer lock infusion connector, color-coded clamps	CSY09	20
Multiple perfusion set for perfusion of three vein grafts with male luer locks, female luer lock infusion connector, color-coded clamps	CSY10	20
Multiple perfusion set for perfusion of three vein grafts with male luer locks and attached vessel cannula, male luer lock infusion connector, color-coded clamps	CSY11	20
Cardioplegia delivery set for retrograde and aortic root perfusion (with vent), female luer lock infusion connector, color-coded clamps	CSY12	20

VENT CATHETERS



Silicone Vent

catalog #	size french	length cm	qty / master carton
SV20	20Fr	35	20

VENT PLUGS



Vent Plug

catalog #	size inch	qty / master carton
VP14	1/4	20
VP38	3/8	20

SAPHENOUS VEIN CANNULAE



Saphenous Vein Cannulae

	catalog #	size french	qty / master carton
vessel cannulae	VC3	1	40
vessel with closed back check valve	VC3B	1	40

SUCTION CONNECTING TUBES



Suction Connecting Tubes

	catalog #	size mm	length mm	qty / master carton
double	TUB69250	6x9	250	40
female funnel	TUB69300	6x9	300	40
	TUB69350	6x9	350	40
	TUB69400	6x9	400	40

ASPIRATION TUBES



catalog #	length cm	qty / master carton
ASP180	180	10
ASP380	380	10


YANKAUER



Yankauer

connecting tubing	catalog #	qty / master carton
bulb tip, non vented	ST01	20
bulb tip, vented	ST02	20
open tip, non vented	ST03	20
open tip, vented	ST04	20

MAMMARY RETRACTOR



description	catalog #
mammary retractor	MR01

The AndoCor mammary retractor is used for the preparation of the internal mammary artery and can be utilised for patients weighing over 25 kg (55lbs).

The mammary retractor is made of surgical stainless steel and has an asymmetrical opening up to 70 mm maximum.

GAS DIFFUSER



catalog #	stylet cm	length cm	qty / master carton
GAS225	40	225	20

SUCKER



rigid sucker

catalog #	size french	length cm	qty / master carton
RS01		33	20



flexible sucker

FS300	17	35	20
-------	----	----	----

HEMOCONCENTRATORS



	description	catalog #	qty / master carton
pre-rinsed	Hemoconcentrator HC250	03HP250	10
	Hemoconcentrator HC600	03HP0600	10
	Hemoconcentrator HC1000	03HP1000	10
	Hemofilter HF2000	03HP2000	10
non-rinsed	Hemoconcentrator HC250	03HC250	10
	Hemoconcentrator HC600	03HC600	10
	Hemoconcentrator HC1000	03HC1000	10
	Hemofilter 2000	03HF2000	10
sets	03HP250, 03HP0600, 03HP1000, 03HC250, 03HC600 and 03HC1000	03HCS001	10
	03HP2000 and 03HF2000	03HFS001	10

These items can also be supplied pre-connected on demand.

HEMOCONCENTRATORS



Technical features

	HC250	HC600	HC1000	HF2000	
Intended use	Paediatric	Adult	Adult	Adult	
Hollow fiber material/Cut off	MediSulfone Polysulfone - 50,000 Daltons	MediSulfone Polysulfone - 50,000 Daltons	MediSulfone Polysulfone - 50,000 Daltons	MediSulfone Polysulfone - 50,000 Daltons	
Potting compound	<i>Polyurethane</i>				
Cartridge housing material	<i>Polycarbonate (PC)</i>				
Cartridge blood ports material	<i>Polycarbonate (PC)</i>				
Cartridge blood connectors	<i>Twist Lock</i>				
Cartridge filtrate/dialysate connectors	Female luer lock	Female luer lock	Female luer lock	Hansen	
Membrane area (m²)	0.25	0.6	1.0	2.0	
Fiber internal diameter (µm)	250	250	250	250	
Fiber external diameter (µm)	350	350	350	350	
Fiber wall thickness (µm)	50	50	50	50	
Fiber length (mm)	140	140	225	280	
Number of fibers	2,400	6,100	6,100	9,800	
Filter housing internal diameter (mm)	24	38	38	49	
Priming blood volume (ml)	18	51	76	146	
Max. transmembrane pressure (mmHg)	500	500	500	500	
Unit length (mm)	180	180	265	325	
Ultrafiltration rate					
<i>Average K_{uf} in vivo 22 ml/h/mmHg/m²</i>	5.5	13.2	22.0	44.0	
Catalog #					
	Pre-rinsed	03HP0250	03HP0600	03HP1000	03HP2000
	Non-rinsed	03HC0250	03HC0600	03HC1000	03HC2000

HEMOFILTER



catalog #	description	qty / master carton
03HF3	adult hemofilter	20
03HF6	adult hemofilter	20
03HF8	adult hemofilter	20

TECHNICAL FEATURES



Catalog #	09HF3	09HF6	09HF8
Ultrafiltration Coefficient	52 ml/hr/mmHg	67 ml/hr/mmHg	78 ml/hr/mmHg
Surface area	1,3m ²	1,6m ²	1,8m ²
Priming volume	80ml	99ml	112ml
Max TMP	500mmHg	500mmHg	500mmHg
Sieving Coefficient Vit B 12	1	1	1
Inulin	1	1	1
B2 – Microglobulin	0,74	0,74	0,74
Albumin	0,01	0,01	0,01
Housing	polycarbonate	polycarbonate	polycarbonate
Potting	polyurethane	polyurethane	polyurethane
Membran type	polyethersulfone	polyethersulfone	polyethersulfone
Fibre Wall thickness	30 µm	30 µm	30 µm
Internal fibre Diameter	200 µm	200 µm	200 µm
Sterilisation	Steam	Steam	Steam

CONNECTORS



Straight Connector

	catalog #	size inch	qty / master carton
without luer lock	04CS001	1/2 x 1/2	20
	04CS002	3/8 x 3/8	20
	04CS003	1/4 x 1/4	20
	04CS004	3/8 x 1/2	20
	04CS005	1/4 x 3/8	20
	04CS006	3/16 x 1/4	20
with luer connectors	04CS0012	1/2 x 1/2	20
	04CS0022	3/8 x 3/8	20
	04CS0032	1/4 x 1/4	20
	04CS0042	3/8 x 1/2	20
	04CS0052	1/4 x 3/8	20
Y-Connector			
	catalog #	size inch	qty / master carton
without luer lock	04CY001	3/8 x 3/8 x 3/8	20
	04CY002	1/4 x 1/4 x 1/4	20
	04CY003	3/8 x 3/8 x 1/2	20
	04CY004	1/4 x 1/4 x 1/4	20
with luer lock	04CY0012	3/8 x 3/8 x 3/8	20

NOTES



Blank area for notes, consisting of several horizontal lines for writing.

ANDOCOR n.v.

Kwikaard 104
B-2980 Zoersel
BELGIUM
T +32 (0)3 380 15 95
F +32 (0)3 380 15 90
info@andocor.com
www.andocor.com