


ANDO COR



# CONTENT

page

## **CARDIOVASCULAR CANNULATION PRODUCTS**

**3 - 11**

---

Arterial cannulae	3
Aortic catheter	7
Venous cannulae	8
Flex line venous cannulae	9
Flex line two stage venous cannulae	10

## **CARDIOPLEGIA PRODUCTS**

**12 - 18**

---

Aortic root cannulae	12
Retrograde cardioplegia cannula	13
Cardioplegia needle	14
Cardioplegia adapters	15
Ostial perfusion cannulae	18

## **VENT CATHETERS**

**19 - 22**

---

Silicone vent catheter	19
PVC vent catheters	20

## **SUCKERS**

**23 - 26**

---

Rigid suckers	23
Flexible suckers	24
Pericardial sump	25
Yankauer suction tube	26

## **GAS DIFFUSER**

**27**

## **HEMOCONCENTRATORS**

**28 - 29**

## **ACCESSORIES**

**30 - 37**

---

Suction connecting tubes	30
Aspiration tubes	30
Porous vent plugs	31
Vented connector cap	31
Vessel cannulae	32
Connectors	33
Cannulation tourniquet sets	34
Quick prime line	35
Tubing organizer	36
Pressure monitoring line	37

## **MAMMARIA RETRACTOR**

**38**

# CARDIOVASCULAR CANNULATION PRODUCTS




## ARTERIAL CANNULAE

Andocor arterial cannulae are intended for use in perfusion of the ascending aorta during cardiopulmonary bypass.



### Features & Benefits

- ✓ Kink resistant reinforced tubing
- ✓ Easy tip orientation thanks to indexing marker
- ✓ Smooth plastic tip



**REINFORCED,  
CURVED TIP WITH FLANGE**

	<i>without vent plug vented red cap</i>	<i>with vent plug</i>	<i>with vented red cap</i>	<i>size french</i>	<i>connection site</i>	<i>length cm</i>	<i>qty / master carton</i>
without connector	A2016	A2016V		20	3/8"	33	20
	B2016	B2016V		20	3/8"	24	20
	A2216	A2216V		22	3/8"	33	20
	B2216	B2216V		22	3/8"	24	20
	A2416	A2416V		24	3/8"	33	20
	B2416	B2416V		24	3/8"	24	20
with connector	A20161		A20161V	20	3/8"	33	20
	B20161		B20161V	20	3/8"	24	20
	A22161		A22161V	22	3/8"	33	20
	B22161		B22161V	22	3/8"	24	20
	A24161		A24161V	24	3/8"	33	20
	B24161		B24161V	24	3/8"	24	20
with luer connector	A20162		A20162V	20	3/8"	33	20
	B20162		B20162V	20	3/8"	24	20
	A22162		A22162V	22	3/8"	33	20
	B22162		B22162V	22	3/8"	24	20
	A24162		A24162V	24	3/8"	33	20
	B24162		B24162V	24	3/8"	24	20

# CARDIOVASCULAR CANNULATION PRODUCTS




## ARTERIAL CANNULAE

Andocor arterial cannulae are intended for use in perfusion of the ascending aorta during cardiopulmonary bypass.



### Features & Benefits

- ✓ Kink resistant reinforced tubing
- ✓ Easy tip orientation thanks to indexing marker
- ✓ Smooth plastic tip



**REINFORCED,  
STRAIGHT TIP WITH RING**

	<i>without vent plug vented red cap</i>	<i>with vent plug</i>	<i>with vented red cap</i>	<i>size french</i>	<i>connection site</i>	<i>length cm</i>	<i>qty / master carton</i>
without connector	A2013	A2013V		20	3/8"	33	20
	B2013	B2013V		20	3/8"	24	20
	A2213	A2213V		22	3/8"	33	20
	B2213	B2213V		22	3/8"	24	20
	A2413	A2413V		24	3/8"	33	20
	B2413	B2413V		24	3/8"	24	20
with connector	A20131		A20131V	20	3/8"	33	20
	B20131		B20131V	20	3/8"	24	20
	A22131		A22131V	22	3/8"	33	20
	B22131		B22131V	22	3/8"	24	20
	A24131		A24131V	24	3/8"	33	20
	B24131		B24131V	24	3/8"	24	20
with luer connector	A20132		A20132V	20	3/8"	33	20
	B20132		B20132V	20	3/8"	24	20
	A22132		A22132V	22	3/8"	33	20
	B22132		B22132V	22	3/8"	24	20
	A24132		A24132V	24	3/8"	33	20
	B24132		B24132V	24	3/8"	24	20

# CARDIOVASCULAR CANNULATION PRODUCTS




## ARTERIAL CANNULAE

Andocor arterial cannulae are intended for use in perfusion of the ascending aorta during cardiopulmonary bypass.



### Features & Benefits

- ✓ Easy tip orientation thanks to indexing marker
- ✓ Smooth plastic tip
- ✓ Kink resistant PVC tubing



**NON-REINFORCED,  
CURVED TIP WITH FLANGE**

	<i>without vent plug vented red cap</i>	<i>with vented red cap</i>	<i>size french</i>	<i>connection site</i>	<i>length cm</i>	<i>qty / master carton</i>
without connector	A2006		20	3/8"	25	20
	A2206		22	3/8"	25	20
	A2406		24	3/8"	25	20
with connector	A20061	A20061V	20	3/8"	25	20
	A22061	A22061V	22	3/8"	25	20
	A24061	A24061V	24	3/8"	25	20
with luer connector	A20062	A20062V	20	3/8"	25	20
	A22062	A22062V	22	3/8"	25	20
	A24062	A24062V	24	3/8"	25	20

# CARDIOVASCULAR CANNULATION PRODUCTS





## ARTERIAL CANNULAE

Andocor arterial cannulae are intended for use in perfusion of the ascending aorta during cardiopulmonary bypass.



### Features & Benefits

- ✓ Easy tip orientation thanks to indexing marker
- ✓ Smooth plastic tip
- ✓ Kink resistant pvc tubing



**NON-REINFORCED,  
STRAIGHT TIP WITH RING**

	<i>without vent plug vented red cap</i>	<i>with vented red cap</i>	<i>size french</i>	<i>connection site</i>	<i>length cm</i>	<i>qty / master carton</i>
without connector	A2003		20	3/8"	25	20
	A2203		22	3/8"	25	20
	A2403		24	3/8"	25	20
with connector	A20031	A20031V	20	3/8"	25	20
	A22031	A22031V	22	3/8"	25	20
	A24031	A24031V	24	3/8"	25	20
with luer connector	A20032	A20032V	20	3/8"	25	20
	A22032	A22032V	22	3/8"	25	20
	A24032	A24032V	24	3/8"	25	20

# CARDIOVASCULAR CANNULATION PRODUCTS




## AORTIC CATHETER

Andocor aortic catheters are intended for use in perfusion of the ascending aorta during cardiopulmonary bypass.



### Features & Benefits

- ✓ Easy tip orientation thanks to indexing marker
- ✓ Movable suture ring
- ✓ Smooth beveled tip



**NON-REINFORCED,  
STRAIGHT TIP WITH RING**

	<i>without vented red cap</i>	<i>with vented red cap</i>	<i>size french</i>	<i>connection site</i>	<i>length cm</i>	<i>qty / master carton</i>
without connector	AC2013		20	3/8"	24	20
	AC2213		22	3/8"	24	20
	AC2413		24	3/8"	24	20
	AC2613		26	3/8"	24	20
	AC3013		30	3/8"	24	20
with connector	AC201341	AC201341V	20	3/8"	24	20
	AC221341	AC221341V	22	3/8"	24	20
	AC241341	AC241341V	24	3/8"	24	20
	AC261341	AC261341V	26	3/8"	24	20
	AC301341	AC301341V	30	3/8"	24	20
with luer connector	AC221342	AC221342V	22	3/8"	24	20
	AC241342	AC241342V	24	3/8"	24	20
	AC261342	AC261342V	26	3/8"	24	20
	AC301342	AC301342V	30	3/8"	24	20

# CARDIOVASCULAR CANNULATION PRODUCTS





## VENOUS CANNULAE

Andocor venous cannulae are intended for drainage from the superior and inferior vena cava.



### Features & Benefits

- ✓ Large selection of sizes
- ✓ Kink resistant reinforced tubing



### **STRAIGHT, REINFORCED WITH BULLET TIP**

<i>straight tip</i>	<i>size french</i>	<i>connection site</i>	<i>length cm</i>	<i>qty / master carton</i>
V2018	20	3/8"	39	20
V2618	26	3/8"	39	20
V3218	32	3/8"	39	20
V3418	34	3/8"	39	20
V3618	36	3/8"	39	20
V4018	40	1/2"	39	20

### **RIGHT ANGLED, REINFORCED WITH BULLET TIP**

<i>right angled tip 90°</i>	<i>size french</i>	<i>connection site</i>	<i>length cm</i>	<i>qty / master carton</i>
V2019	20	3/8"	39	20
V2619	26	3/8"	39	20
V3219	32	3/8"	39	20
V3419	34	3/8"	39	20
V3619	36	3/8"	39	20
V4019	40	1/2"	39	20

### **NON-REINFORCED**

<i>straight tip</i>	<i>size french</i>	<i>connection site</i>	<i>length cm</i>	<i>qty / master carton</i>
VC3208	32	3/8"	40	20



# CARDIOVASCULAR CANNULATION PRODUCTS



## FLEX LINE VENOUS CANNULAE

Andocor venous cannulae are intended to drain venous blood from the right atrium and caval veins.



### Features & Benefits

- ✓ Open lighthouse tip
- ✓ Single piece construction

### **STRAIGHT, REINFORCED WITH LIGHTHOUSE TIP**

<i>straight tip</i>	<i>size french</i>	<i>connection site</i>	<i>length cm</i>	<i>qty / master carton</i>
01V201L8	20	3/8"	39	20
01V221L8	22	3/8"	39	20
01V241L8	24	3/8"	39	20
01V261L8	26	3/8"	39	20
01V281L8	28	3/8"	39	20
01V301L8	30	3/8"	39	20
01V321L8	32	3/8"	39	20
01V341L8	34	3/8"	39	20
01V361L8	36	1/2"	39	20

### **RIGHT ANGLED, REINFORCED WITH LIGHTHOUSE TIP**

<i>right angled tip 90°</i>	<i>size french</i>	<i>connection site</i>	<i>length cm</i>	<i>qty / master carton</i>
01V201L9	20	3/8"	39	20
01V221L9	22	3/8"	39	20
01V241L9	24	3/8"	39	20
01V261L9	26	3/8"	39	20
01V281L9	28	3/8"	39	20
01V301L9	30	3/8"	39	20
01V321L9	32	3/8"	39	20
01V341L9	34	3/8"	39	20
01V361L9	36	1/2"	39	20

# CARDIOVASCULAR CANNULATION PRODUCTS



## FLEX LINE TWO STAGE VENOUS CANNULAE

Andocor two stage venous cannulae are intended for use in venous drainage via the right atrium and inferior vena cava simultaneously during cardiopulmonary bypass surgery up to six hours or less.



### Features & Benefits

- ✓ Open and flexible lighthouse tip
- ✓ Single piece construction
- ✓ Easy positioning
- ✓ Obturator
- ✓ Insertion depth marker



**PROXIMAL AND DISTAL  
REINFORCED  
WITH LIGHTHOUSE TIP**

	<i>standard body</i>	<i>flat body</i>	<i>size french</i>	<i>length cm</i>	<i>qty / master carton</i>
without connector	01V32L40	01V32L40FB	32/40	40	20
	01V34L46	01V34L46FB	34/46	40	20
	01V36L46	01V36L46FB	36/46	40	20
	01V36L51	01V36L51FB	36/51	40	20
with connector	01V32L401	01V32L401FB	32/40	40	20
	01V34L461	01V34L461FB	34/46	40	20
	01V36L461	01V36L461FB	36/46	40	20
	01V36L511	01V36L511FB	36/51	40	20
with luer lock	01V32L402	01V32L402FB	32/40	40	20
	01V34L462	01V34L462FB	34/46	40	20
	01V36L462	01V36L462FB	36/46	40	20
	01V36L512	01V36L512FB	36/51	40	20

# CARDIOVASCULAR CANNULATION PRODUCTS




## FLEX LINE TWO STAGE VENOUS CANNULAE

Andocor two stage venous cannulae are intended for use in venous drainage via the right atrium and inferior vena cava simultaneously during cardiopulmonary bypass surgery up to six hours or less.



### Features & Benefits

- ✓ Bullet tip
- ✓ Single piece construction
- ✓ Easy positioning
- ✓ Obturator
- ✓ Insertion depth marker



**PROXIMAL AND DISTAL  
REINFORCED  
WITH BULLET TIP**



	<i>standard body</i>	<i>flat body</i>	<i>size french</i>	<i>length cm</i>	<i>qty / master carton</i>
without connector	01V34146	01V34146FB	34/46	40	20
	01V36146	01V36146FB	36/46	40	20
	01V36151	01V36151FB	36/51	40	20
with connector	01V341461	01V341461FB	34/46	40	20
	01V361461	01V361461FB	36/46	40	20
	01V361511	01V361511FB	36/51	40	20
with luer lock	01V341462	01V341462FB	34/46	40	20
	01V361462	01V361462FB	36/46	40	20
	01V361512	01V361512FB	36/51	40	20

# CARDIOPLEGIA PRODUCTS



## AORTIC ROOT CANNULAE

Andocor aortic root cannulae are intended for use during cardiopulmonary bypass for antegrade delivery of cardioplegia solutions and venting of the heart. The cannula may also be used to aspirate air from the aorta at the conclusion of the cardiac procedure.



### + Features & Benefits

- ✓ With & without vent line
- ✓ Different configurations

	<i>code</i>	<i>size needle</i>	<i>size french</i>	<i>length cm</i>	<i>qty / master carton</i>
without vent line	ARC07	18G	7	15	20
	ARC04	16G	8	15	20
	ARC01	14G	10	15	20
	ARC14	16G	8	28	20
with vent line	ARC08	18G	7	15	20
	ARC05	16G	8	15	20
	ARC02	14G	10	15	20

# CARDIOPLEGIA PRODUCTS




## RETROGRADE CARDIOPLEGIA CANNULA

Andocor retrograde cardioplegia cannulae are intended to be connected to the cardioplegia line in order to infuse cardioplegic solution and blood into the patient's coronary sinus during open heart procedure.



### Features & Benefits

- ✓ Self-inflating
- ✓ Silicone balloon
- ✓ Guide wire stylet



<i>code</i>	<i>balloon diameter</i>	<i>size french</i>	<i>length cm</i>	<i>qty / master carton</i>
RCC01	18 mm	14	31	20

# CARDIOPLEGIA PRODUCTS





## CARDIOPLEGIA NEEDLE

Andocor cardioplegia needles are intended for use in delivery of cardioplegia solution during the cardiopulmonary bypass surgery.



### Features & Benefits

- ✓ Protective cap
- ✓ Stainless steel needle
- ✓ Luer lock connection



	<i>code</i>	<i>total length cm</i>	<i>needle length cm</i>	<i>size needle</i>	<i>qty/ master carton</i>
standard	CN16	5	1,4	16G	20

# CARDIOPLEGIA PRODUCTS




## CARDIOPLEGIA ADAPTERS

Andocor cardioplegia adaptors connect to the ARC or vessel cannulae and are intended for the delivery of cardioplegia solution or venting the heart during cardiopulmonary bypass.



### Features & Benefits

- ✓ Different configurations
- ✓ Color-coded clamps



CSY03 perfusion adapter with 2 male luer lock and 1 female luer lock and color-coded clamps

CSY03



CSY04 multiple perfusion adapter for aortic root and three vein grafts

CSY04



CSY05 multiple perfusion adapter for aortic root and three vein grafts with attached vessel cannulae

CSY05




CSY06 dual perfusion set with female luer lock, cone connector and color coded clamps

CSY06











CSY07 multiple perfusion adapter with vent line and vessel cannula

CSY07



# CARDIOPLEGIA PRODUCTS











<p>CSY08</p>	<p>aortic perfusion set with vent (long = 20cm)</p>	<p>CSY08</p> 
<p>CSY09</p>	<p>dual perfusion set with female luer lock, male luer lock connector, color coded clamps</p>	<p>CSY09</p> 
<p>CSY10</p>	<p>multiple set for perfusion of three vein grafts</p>	<p>CSY10</p> 
<p>CSY11</p>	<p>multiple set for perfusion of three vein grafts with attached vessel cannulae</p>	<p>CSY11</p> 
<p>CSY12</p>	<p>multiple perfusion adapter for retrograde and aortic root with vent</p>	<p>CSY12</p> 
<p>CSY13</p>	<p>aortic perfusion set with vent (short = 11cm)</p>	<p>CSY13</p> 
<p>CSY14</p>	<p>aortic perfusion set with 1/4x1/4x1/4 connector</p>	<p>CSY14</p> 
<p>CSY15</p>	<p>straight adapter, 14cm</p>	<p>CSY15</p> 



# CARDIOPLEGIA PRODUCTS



CSY16	straight adapter, 15cm	
CSY17	perfusion adapter 1/4"	
CSY18	perfusion adapter 3/8"	
CSY19	perfusion adapter 1/4" with 3-way stopcock	
CSY20	perfusion adapter with 2 female luer lock and 1 male luer lock and color-coded clamps	
CSY21	recirculation y-adapter	
CSY22	multiple perfusion adapter with vent line	
CSY23	six multiple perfusion adapter	

# CARDIOPLEGIA PRODUCTS





## OSTIAL PERFUSION CANNULAE

Andocor ostial perfusion cannulae are intended for use in delivery of cardioplegia solution directly to the coronary arteries during cardiopulmonary bypass surgery.



### Features & Benefits

- ✓ Basket tip
- ✓ Different configurations
- ✓ Malleable stainless steel shaft



	<i>code</i>	<i>size french</i>	<i>length cm</i>	<i>qty / master carton</i>
with 45° tip	OC1045	10	21	20
	OC1245	12	21	20
	OC1445	14	21	20
with 90° tip	OC1090	10	21	20
	OC1290	12	21	20
	OC1490	14	21	20

# VENT CATHETER




## SILICONE VENT CATHETER

Andocor vent catheters are medical devices to decompress or vent the heart during cardiopulmonary bypass cannulation. Vent catheters are used as a suction catheter to assist in the complete evacuation of blood from the heart.



### Features & Benefits

- ✓ With guide wire stylet
- ✓ Soft and flexible silicone tubing



**VENTRICULAR VENT  
CATHETER, SILICONE,  
WITH CLOSED TIP**

<i>code</i>	<i>size french</i>	<i>tip perforation cm</i>	<i>length cm</i>	<i>qty / carton</i>
SV20	20	6,5	35	20

# VENT CATHETER




## PVC VENT CATHETER

Andocor vent catheters are medical devices to decompress or vent the heart during cardiopulmonary bypass cannulation. Vent catheters are used as a suction catheter to assist in the complete evacuation of blood from the heart.



### Features & Benefits

- ✓ Variety of configurations
- ✓ Open tip



**VENTRICULAR VENT  
CATHETER, PVC,  
WITH OPEN TIP, STRAIGHT**

	<i>code</i>	<i>size french</i>	<i>tip perforation cm</i>	<i>length cm</i>	<i>qty / carton</i>
without stylet, 1/4" connector	PVSO171	17	5,5	32	20
with removable stylet, 1/4" connector	PVSOR181	18	6	32	20
with internal stylet, 1/4" connector	PVSOI181	18	6	32	20
with guide wire stylet, 1/4" connector	PVSOG181	18	6	32	20

# VENT CATHETER





## PVC VENT CATHETER

Andocor vent catheters are medical devices to decompress or vent the heart during cardiopulmonary bypass cannulation. Vent catheters are used as a suction catheter to assist in the complete evacuation of blood from the heart.



### Features & Benefits

- ✓ Variety of configurations
- ✓ Closed tip



**VENTRICULAR VENT  
CATHETER, PVC,  
WITH CLOSED TIP, STRAIGHT**

	<i>code</i>	<i>size french</i>	<i>tip perforation cm</i>	<i>length cm</i>	<i>qty / carton</i>
with stylet, 1/4" connector	PVSB181	18	6	32	20
with guide wire stylet, 1/4" connector	PVSBG181	18	6	32	20
with removable stylet, 1/4" connector	PVSBR181	18	6	32	20
with internal stylet, 1/4" connector	PVSB181	18	6	32	20
with guide wire stylet, 1/4" connector	PVSBG201	20	6	32	20
with removable stylet, 1/4" connector	PVSBR201	20	6	32	20

# VENT CATHETER





## PVC VENT CATHETER

Andocor vent catheters are medical devices to decompress or vent the heart during cardiopulmonary bypass cannulation. Vent catheters are used as a suction catheter to assist in the complete evacuation of blood from the heart.



### Features & Benefits

- ✓ Closed tip
- ✓ Depth marker
- ✓ Different sizes



**VENTRICULAR VENT  
CATHETER, PVC, WITH  
CLOSED TIP, RIGHT ANGLED**

	<i>code</i>	<i>size french</i>	<i>tip perforation cm</i>	<i>length cm</i>	<i>qty / carton</i>
without stylet	PVBA13	13	3,8	32	20
	PVBA20	20	3,8	36,5	20

# SUCKERS





## RIGID SUCKERS

Andocor suckers are used for the removal of excess fluid from the surgical field.



### Features & Benefits

- ✓ Metal suction tube with plastic handle
- ✓ Easy & comfortable directing



	<i>code</i>	<i>length cm</i>	<i>qty / carton</i>
Rigid tip	RS01	33	20

# SUCKERS





## FLEXIBLE SUCKERS

Andocor suckers are used for the removal of excess fluid from the surgical field.



### Features & Benefits

- ✓ Provide continuous suction
- ✓ Flexible or malleable tubing



**FLEXIBLE / MALLEABLE  
INTRACARDIAC SUCKERS  
WITH FLUTED TIP**

	<i>code</i>	<i>size french</i>	<i>length cm</i>	<i>connection site</i>	<i>qty / carton</i>
flexible tubing	FS300	17	36	1/4"	20
malleable tubing	FS301	17	45	1/4"	20



# SUCKERS




## PERICARDIAL SUMP

Andocor suckers are designed to remove fluid from the surgical field.



### Features & Benefits

- ✓ Wire-wound tip
- ✓ Drainage of pericardial sack



**REINFORCED TIP**

	<i>code</i>	<i>size french</i>	<i>length cm</i>	<i>connection site</i>	<i>qty / carton</i>
Reinforced tip	PS40	16	36	1/4"	20



### Features & Benefits

- ✓ Smooth and flush design
- ✓ Metal tip
- ✓ External diameter of 6 mm



**METAL TIP**

	<i>code</i>	<i>size french</i>	<i>length cm</i>	<i>connection site</i>	<i>qty / carton</i>
Metal tip	PS40M	18	36	1/4"	20

# SUCKERS





## YANKAUER SUCTION TUBE

Andocor yankauer suction tubes are intended to aspirate body waste fluid during surgery.



### Features & Benefits

- ✓ Different configurations



	<i>tip</i>	<i>code</i>	<i>qty / master carton</i>
vented	open	ST04	20
	bulb	ST02	20
non-vented	open	ST03	20
	bulb	ST01	20

# GAS DIFFUSER




## GAS DIFFUSER

Andocor gas diffuser is a medical device for intraoperative wound ventilation with CO<sub>2</sub> to prevent air embolism and bacteriological contamination during surgery.



### Features & Benefits

- ✓ Innovative microporous diffuser tip
- ✓ Bacterial filter
- ✓ Hydrophobic diffuser tip remains functional even when wet
- ✓ Malleable section allows precise placement



<i>code</i>	<i>stylet cm</i>	<i>length cm</i>	<i>qty / master carton</i>
GAS225	40	225	20

# HEMOCONCENTRATORS



## HEMOCONCENTRATORS

Andocor hemoconcentration filters are intended to relieve or mitigate overhydration in patients undergoing cardiopulmonary procedures and to increase concentrations of cells and proteins in the blood.



### Features & Benefits

- ✓ Pre-rinsed or non-rinsed filters
- ✓ Standard or custom-made blood lines



	<i>description</i>	<i>code</i>	<i>qty / master</i>
<i>pre-rinsed</i>	Hemoconcentrator HC0250	03HP0250	10
	Hemoconcentrator HC0600	03HP0600	6
	Hemoconcentrator HC1000	03HP1000	6
	Hemofilter HF2000	03HP2000	6
<i>non-rinsed</i>	Hemoconcentrator HC0250	03HC0250	10
	Hemoconcentrator HC0600	03HC0600	6
	Hemoconcentrator HC1000	03HC1000	6
	Hemofilter 2000	03HF2000	6
<i>sets</i>	03HP0250, 03HP0600, 03HP1000, 03HC250, 03HC0600 and 03HC1000	03HCS001	40
	03HP2000 and 03HF2000	03HFS002	40

# HEMOCONCENTRATORS



## HEMOCONCENTRATORS

Andocor hemoconcentration filters are intended to relieve or mitigate overhydration in patients undergoing cardiopulmonary procedures and to increase concentrations of cells and proteins in the blood.

## TECHNICAL FEATURES

	<b>HC250</b>	<b>HC600</b>	<b>HC1000</b>	<b>HF2000</b>
<i>Intended use</i>	<i>Paediatric</i>	<i>Adult</i>	<i>Adult</i>	<i>Adult</i>
<i>Hollow fiber material/Cut off</i>	<i>MediSulfone Polysulfone - 50,000 Daltons</i>	<i>MediSulfone Polysulfone - 50,000 Daltons</i>	<i>MediSulfone Polysulfone - 50,000 Daltons</i>	<i>MediSulfone Polysulfone - 50,000 Daltons</i>
<i>Potting compound</i>		Polyurethane		
<i>Cartridge housing material</i>		Polycarbonate (PC)		
<i>Cartridge blood ports material</i>		Polycarbonate (PC)		
<i>Cartridge blood connectors</i>		Twist Lock		
<i>Cartridge filtrate/dialysate connectors</i>	Female luer lock	Female luer lock	Female luer lock	Hansen
<i>Membrane area (m<sup>2</sup>)</i>	0.25	0.6	1.0	2.0
<i>Fiber internal diameter (µm)</i>	250	250	250	250
<i>Fiber external diameter (µm)</i>	350	350	350	350
<i>Fiber wall thickness (µm)</i>	50	50	50	50
<i>Fiber length (mm)</i>	140	140	225	280
<i>Number of fibers</i>	2,4	6,1	6,1	9,8
<i>Filter housing internal diameter (mm)</i>	24	38	38	49
<i>Priming blood volume (ml)</i>	18	51	76	146
<i>Max. transmembrane pressure (mmHg)</i>	500	500	500	500
<i>Unit length (mm)</i>	180	180	265	325
<i>Ultrafiltration rate</i>				
<i>Average Kuf in vivo 22 ml/h/mmHg/m<sup>2</sup></i>	5.5	13.2	22.0	44.0
<i>Pre-rinsed</i>	03HP0250	03HP0600	03HP1000	03HP2000
<i>Non-rinsed</i>	03HC0250	03HC0600	03HC1000	03HF2000

# ACCESSORIES




## SUCTION CONNECTING TUBES

Andocor suction connecting tubes are used for removal of excess fluid from a surgical field.



### Features & Benefits


- ✓ Double female funnel
- ✓ Double packed



<i>code</i>	<i>size mm</i>	<i>length mm</i>	<i>qty / master carton</i>
TUB69180	6x9	180	20
TUB69380	6x9	380	10

## ASPIRATION TUBES

Andocor aspiration tubes are used as an accessory for blood aspiration during cardiovascular surgery and cannulation.




<i>code</i>	<i>length cm</i>	<i>qty / master carton</i>
ASP180	180	25
ASP380	380	25

# ACCESSORIES




## POROUS VENT PLUGS



<i>code</i>	<i>size</i>	<i>qty / master carton</i>
VP38	3/8"	40

## VENTED CONNECTOR CAP



<i>code</i>	<i>size</i>	<i>qty / master carton</i>
VCC01	3/8"	40

# ACCESSORIES





## VESSEL CANNULAE

Andocor vessel cannulae facilitate administration of cardioplegia solution to the heart.



### Features & Benefits

- ✓ Choice of flexible or rigid tube



	<i>code</i>	<i>length cm</i>	<i>qty / master carton</i>
vessel cannulae, transparant, flexible	VC3	5	40
vessel cannulae with closed back check valve	VC3B	6,5	40
vessel cannulae, white, rigid	VC04	5	40



# ACCESSORIES




## CONNECTORS

Andocor connectors are accessories for the connection of cannulae to the extracorporeal cardiopulmonary bypass system during cardiovascular surgery and cannulation.



### Features & Benefits

- ✓ Wide range



	code	outlet 1	size inch outlet 2	outlet 3	qty / master carton	
<b>straight connector</b> <i>without luer lock</i>	04CS001	1/2"	1/2"		20	
	04CS002	3/8"	3/8"		20	
	04CS003	1/4"	1/4"		20	
	04CS004	3/8"	1/2"		20	
	04CS005	1/4"	3/8"		20	
	04CS006	3/16"	1/4"		20	
	04CS007	3/16"	3/16"		20	
	04CS008	3/8"	3/16"		20	
	<i>with luer connectors</i>	04CS0012	1/2"	1/2"		20
		04CS0022	3/8"	3/8"		20
		04CS0032	1/4"	1/4"		20
		04CS0042	3/8"	1/2"		20
		04CS0052	1/4"	3/8"		20
		04CS0062	3/16"	1/4"		20
		04CS0072	3/16"	3/16"		20
<b>Y-connector</b> <i>without luer lock</i>		04CY001	3/8"	3/8"	3/8"	20
	04CY002	1/4"	1/4"	1/4"	20	
	04CY003	3/8"	3/8"	1/2"	20	
	04CY004	1/4"	1/4"	3/8"	20	
	04CY005	1/2"	1/2"	1/2"	20	
	04CY006	3/8"	1/4"	3/8"	20	
	04CY007	1/2"	1/4"	1/2"	20	
	<i>with luer lock</i>	04CY0012	3/8"	3/8"	3/8"	20

# ACCESSORIES





## CANNULATION TOURNIQUET SETS

Andocor cannulation tourniquet sets are used to secure arterial and venous cannulae during CBP.



### Features & Benefits

- ✓ Wide range of configurations




<i>code</i>	<i>snare</i>	<i>red</i>	<i>blue</i>	<i>clear</i>	<i>length cm</i>	<i>qty / master carton</i>
CT01	1	1	1	3	17,6	20
CT02	1	1	2	3	17,6	20
CT03	2	1	2	5	17,6	20
CT04	1	1	1		17,6	20
CT05	1	2	1		17,6	20
CT06	2	2	2	4	17,6	20
CT07	2	4	2	2	17,6	20
CT08	1	5			17,6	20
CT09	1			3	17,6	20

# ACCESSORIES



## QUICK PRIME LINE

Andocor quick prime line facilitates the transfer of fluids during the priming of the extracorporeal circuit.




<i>code</i>	<i>length cm</i>	<i>qty / master carton</i>
QPL100	100	20

# ACCESSORIES



## TUBING ORGANIZER



<i>code</i>	<i>slots &amp; type</i>	<i>qty</i>	<i>qty/ master carton</i>
TO10	10 tubes		20
	1/8"	1	
	1/4"	6	
	3/8"	2	
	1/2"	1	

# ACCESSORIES



## PRESSURE MONITORING LINE

Andocor pressure monitoring lines are used for extending medical tubing and pressure monitoring.



<i>code</i>	<i>tube end</i>	<i>length cm</i>	<i>qty / master carton</i>
PML180MM	male/male	180	20
PML180MF	male/female	180	20

# ACCESSORIES




## MAMMARIA RETRACTOR

The mammaia retractor is used for the preparation of the internal mammary artery and can be utilised for patients weighing over 25 kg (55 lbs).



### Features & Benefits

- ✓ Surgical stainless steel
- ✓ Asymmetrical opening up to 70 mm max



*code*

MR01

*qty /  
master carton*

1



**ANDOCOR n.v.**

Kwikaard 104  
B-2980 Zoersel  
BELGIUM  
T +32 (0)3 380 15 95  
F +32 (0)3 380 15 90  
[info@andocor.com](mailto:info@andocor.com)  
[www.andocor.com](http://www.andocor.com)